



# PEDIATRICS

product catalog

Your Logo Here



### The MED Group Gold Seal of Approval

This catalog is filled with MED Preferred products identified by The MED Group Gold Seal of Approval. Made by leading healthcare manufacturers, these products are fully backed by the manufacturer's warranty.

MED Preferred products are chosen based on quality and value, do not settle for anything less. MED Group Members provide professional

service before and after the sale to insure the products purchased meet your needs today as well as in the future. Please call your local MED Group Member and experience a higher level of service and professionalism.

Please utilize this catalog and ask your MED Group Member for assistance in helping you locate the products right for your application.

*created by:*



[www.medgroup.com](http://www.medgroup.com)  
800.825.5633

**All products in this book may not be available from all providers.**

# Table of Contents

4 - 5  
POWER WHEELCHAIRS

11 - 12  
TOILETING  
& BATHING

6 - 7  
MANUAL WHEELCHAIRS

13 - 16  
POSITIONING CHAIRS  
& SEATS

8 - 10  
STROLLERS  
& POSITIONING

17 - 19  
STANDERS



# Power Wheelchairs



## Quantum® Rehab 6000Z

- Terrain busting torque and enhanced speed
- Rugged, sporty styling
- Superb climbing and transitioning capability
- Solid track record of durability and customer satisfaction
- Highly maneuverable
- Compatible with an impressive array of electronics and seating and positioning systems to keep up with a growing child's needs
- Available high speed motor package



## Quantum Rehab Dynamo® 1107 featured with Pediatric Tilt

- The Dynamo 1107 is the perfect power chair choice for pediatric use.
- With its small, compact design that fits easily into kid-sized spaces and a 20.5" turning radius, the Dynamo 1107 allows kids to keep pace with their friends at home, at school and on the playground.
- The Dynamo 1107 features direct-connect components that eliminate wiring harnesses for the fastest and easiest assembly and disassembly, helping to keep kids active and on the go.



# Power Wheelchairs



## Quantum® Rehab Pediatric 600

- Compact, patent pending Mid-Wheel 6 design and powerful GLIDE motors let the Quantum 600 keep up with today's active kids.
- The Quantum 600 is versatile enough to meet even the most unique pediatric requirements while delivering the high performance needed to keep up with a jam-packed schedule.



## Quantum® Rehab Pediatric 610

- The Quantum 610's highly compact size provides easy access to kid-sized spaces combined with the true all-around performance of our patent pending Mid-Wheel 6 design.
- GLIDE motors and excellent rehab capability make the Quantum 610 a powerful and complete pediatric power chair option.



# Manual Wheelchairs



## Sunrise Medical Zippie®

- Easy-pushing Zippie wheelchair is both rigid and adaptable.
- Rigid frame (non-folding) construction is durable and strong.
- Open frame design allows for easy integration of seating systems.
- Features a wide range of adjustments and swing-away footrests.
- Designed with growth in mind, each chair comes with one free growth kit.



## Sunrise Medical Zippie® 2

- Lightweight, folding chair designed just for children.
- Locking brace gives it an easy-pushing, solid feel.
- Extremely adjustable and has the ability to grow with a child.
- Swing-away footrests and quick-release components make it easy to transport.
- One time free Growth Kit allows the chair to be 'grown' free of charge.
- One-Arm Drive Option – Mounts on the right or left side, is easily folded and weighs less than 5 lbs.



## Sunrise Medical Zippie® GS

- Designed to easily 'grow' without the need for new parts, the Zippie GS grows with them.
- Built-in growth system provides the greatest range of width, depth and angle adjustment available today.
- Unique side frames add even more flexibility and stability.
- Proven design and the highest quality construction for trouble-free durability.
- An all new universal axle plate.
- Adjustable back upholstery allows for a contoured fit without adding the weight and complexity of some seating solutions.
- Various choices to meet any mobility need: rigid, folding or reverse configuration.

# Manual Wheelchairs



## Sunrise Medical Zippie® TS

- Built-in frame expansion.
- Both rigid and folding frames.
- Simple system for tilting-in-space up to 45 degrees for proper positioning and comfort.
- Lightweight welded aluminum construction and folding frame option for portability.
- Designed to easily accept the widest possible variety of seating solutions.
- Back post mounting system adjusts for proper seat depth, seat-to-back angle and folds for portability.
- All of the latest in pediatric accessories.



## Sunrise Medical Zippie® IRIS

- Based on the patented rotation-in-space design of the Quickie IRIS, downsized to fit the pediatric user.
- A smaller rotational arc to better match a pediatric user's center of gravity.
- 60 degrees tilt advantage – as the seat frame rotates past 45 degrees, a significant amount of pressure transfers from the seat to the back.
- Features 3 tilt ranges: -5 to 50 degrees, 0 to 55 degrees and 5 to 60 degrees.
- Low seat-to-floor height, standard 13" front STFH, with drop seat, as low as 11 in.
- A transit option is available for 150 lb. user weight-limits.



## Sunrise Quickie® Kidz Bike

- The Quickie® Kidz Bike is the perfect addition to your child's recreational activities, allowing them the freedom of biking outside with all their friends.
- Transferring is easy with the sliding seat and flip-up footrest on Kidz Bike.
- Sliding seat and adjustable depth accommodate for growth and riders of different sizes.
- At a seat height of 17", Children on the Quickie® Kidz Bike are visible to vehicles and pedestrians.
- To accommodate for different disabilities, the back angle is adjustable to fit a variety of needs.
- A number of color choices and a single or three speed gear.

# Strollers



## Sammons Preston Special Tomato® Jogger

- Lightweight, maneuverable stroller that does not compromise positioning.
- Super-light aluminum frame.
- Comfortable padded seat reclines and protective windowed canopy.
- Soft-Touch for children who need minimal to moderate positioning.
- Liners easily attach and folds with the stroller.



## Sammons Preston Maclaren® Major Positioning Push Chair

- Developed especially for special needs children to meet the needs of parents for an everyday portable stroller to transport their child.
- The sturdy aluminum chassis and 7" wheels accommodates children up to 140 lbs.
- Swivel-front wheels for easy maneuvering and fixed back wheels with brakes.
- 3-position, detachable footrests.
- Folds like an umbrella stroller for easy transport.
- Optional Accessories: shopping bag, padded seat, sunshade, lateral inserts, chest pad, wind/rain shield and stroller boot.







## Otto Bock Kimba® Cross

- Best pediatric positioning jogging stroller in the market for health conscious caregivers.
- Low center of gravity.
- Large pneumatic tires.
- Easy handling even in rugged terrain.



## Otto Bock Eco-Buggy™

- Powder-coated aluminum frame resist scratching and rust, holds loads up to 110lbs.
- The Eco-Buggy™ features twin casters and rear wheels with an easily accessible and solid footbrake.
- The front casters come standard with swivel wheel locks for easier maneuvering on rough terrain.
- Standard swing-away footrests make transferring the child a breeze.
- The seat and back upholstery are reinforced to ensure stability. The back angle can be reclined for varying hip angle positions.
- Seat depth can be adjusted from 12" to 14" to accommodate growth.
- Passed WC-19 crash test and is bus transportable.



# Strollers



## Stealth Products FLASH (Economy)

- Durable and lightweight seat and back upholstery.
- Back angle adjustment of 85, 90 and 95 degrees.
- Offset arms provide additional width.
- Seat depth adjustments.
- Single platform front rigging – extension adjustable.
- Folds for convenient transport and storage.



## Stealth Products LIGHTNING (Performance)

- Heavy duty upholstery includes moldable 'stays'.
- Back angle adjustment of 85, 90 and 95 degrees.
- Offset arms provide additional width.
- Seat depth adjustments.
- Storage basket standard.
- Independent swing-away front riggings with angle, lateral and fore/aft adjustments.
- Crash tested – optional Tie Down Kit available.
- Folds for convenient transport and storage.



## Stealth Products SPRITE (Seating & Positioning)

- Solid seats and backs from Stealth Products.
- Accommodates most aftermarket seats and backs.
- Seat angle adjustable to 5, 10 and 20 degrees.
- Back angle adjustment of 85, 90 and 95 degrees.
- Offset arms provide additional width.
- Seat depth adjustments.
- Storage basket standard.
- Independent swing-away front riggings with angle, lateral and fore/aft adjustments.
- Crash-tested – optional Tie Down Kit available.
- Folds for convenient transport and storage.

# Toileting & Bathing



## Sammons Preston Leckey® Advanced Bath Chair

- Reaches new levels of functionality and adjustability inside and outside the tub.
- Features fully adjustable chair legs, articulating back and leg support, lightweight but strong stainless steel frame and folds flat for easy storage or transport.
- All-around padded chest and hip belt supports result in superior positioning comfort.
- Easily adjusts for low and high tub positioning.
- Includes non-slip foam pads for secure tub placement and protection.
- Choose from 4 sizes to fit the youngest child through adolescence.

## Sammons Preston Leckey® Advanced Shower Trolley

- Simplifies transfers from the tub to the shower.
- Easily attaches to the Advanced Bath Chair.
- Brings child to the ideal height within a strong, stable, maneuverable base.



## Sammons Preston Leckey® Basic Bath Support

- Provides the same adjustable positioning angles as the Advanced Shower Chair, but at a fixed height.
- One size.
- Includes hip belt & chest belt.



# Toileting & Bathing



## Sammons Preston Leckey® Toilet/Shower Chair

- 3 sizes fit a wide range of children and adolescents.
- This versatile shower chair also works over the toilet or as a stand-alone commode.
- Among its many features are a reversible seat back for upright or slightly prone positioning, fully padded seat back frame, removable pommel/splash guard and 4 lockable swivel casters.



## Sammons Preston Leckey® Potty Trainer

- This attractive toileting solution does not compromise the positional needs of the less-involved child.
- Incorporates many of the features of the Easy Seat to ensure comfortable, confident seating.
- Features include molded ABS plastic commode bucket with splash guard, height and angle-adjustment backrest, 45 degrees hip belt, adjustable side pads, height-adjustable armrests and independently height-adjustable legs.
- Includes non-wheeled, adjustable legs, hip belt and 1" side pads.
- 4 sizes accommodate small children through adolescents.



# Positioning Chairs & Seats

## Sammons Preston Special Tomato Soft-Touch™ Therapy Rolls & Wedges



- Provides soft comfort and support for a variety of therapeutic positioning options.
- Extremely durable and easy to clean.
- Latex-free, seamless material has antimicrobial agent.
- Resists peels and tears, impermeable to fluids.
- Stand-on-End design simplifies storage.
- 3 year manufacturer's warranty.

## Sammons Preston Special Tomato Soft-Touch™ Sitters



- Supportive sitters offer affordable alternative for positioning at home, school or out in the community.
- Same features as Soft-Touch Therapy Rolls & Wedges.
- Modular cushion seating surface allows for personal adjustment and years of growth.
- Each of the 3 sizes can be securely attached to standard chairs.
- Choose the Sitter only or the Sitter with special options such as the Mobile Floor Kit and the Mobile Tilt and the Wedge Kit shown here.

## Sammons Preston Special Tomato Multi-Positioning Seat



- Choose from 3 systems: Multi-Positioning Seat (MPS), MPS Mobile Floor Kit and MPS Push Chair Kit.
- All include adjustable anti-thrust seat, head and back cushions, seat shell and adjustable footrest.
- The first kit comes with a Mobile Floor Base, and the second has a bus transportable Push Chair Base.
- Each system can be customized with optional accessories such as Standard Headrest, Headrest with Laterals, Mobile Floor Base, Extended Seat Cushion or Tray.
- No matter which system you choose, count on using it everywhere – at home, in the care, at school, in the school bus or anywhere your child needs to go!

# Positioning Chairs & Seats



## Sammons Preston Leckey® Mygo Seating System

- A flexible sacral cushion designed to position the pelvis in a neutral or slightly anterior tilt.
- A pelvic positioning harness with adjustable hip guides helps maintain medial alignment.
- Independently adjustable upper leg supports accommodate leg length discrepancy.
- Optional multi-positional foot plates allow the feet to be positioned behind the knees.
- Anatomically correct backrest pivot point is adjustable in depth and height to accommodate posture and growth.



## Sammons Preston Leckey® Easy Seat Complete

- This extremely attractive seating solution grows with the child without compromising positioning needs.
- Incorporates adjustments for seat depth, leg height and armrest height.
- Includes: tray, 1" side pads, foot plate with straps, wheeled legs and lap strap/hip belt.
- 4 sizes accommodate children through adolescents.



# Positioning Chairs & Seats



## Sammons Preston Special Tomato Liners

- Soft and supportive, providing comfortable and affordable alternative seating.
- Made of a latex-free seamless material that are peel and tear-resistant, impermeable to fluids and soft yet durable.
- Withstands cracking, is antimicrobial and cleans easily with common disinfectants.
- Great for school and home use. Use for positioning in strollers or chairs.



## Sammons Preston Special Tomato Height Rite Chair

- Stylish, comfortable seating at home or school. Natural hardwood encourages interaction with peers, while providing minimal to moderate support for the child.
- Has a solid, adjustable footrest that gives the child a safe and secure platform for proper positioning of the legs and feet.
- Seat and footrest adjust in height and depth, allowing the chair to "grow" with the child throughout the years.



## Sammons Preston Leckey Sleep Form

- Provides a perfect night sleep by supporting a symmetric posture. Forms to the individual contours of children, holding them in a secure and comfortable position in their preferred sleeping posture.
- Can be used in a cot or bed and are designed specifically for children within three different age groups.
- System is comfortable, supportive, flexible, easy to use and can be easily transported.

# Positioning Chairs & Seats



## Sammons Preston Leckey® Squiggles Seating System

- Designed for children aged from 1 to 5 years to provide them with comfortable and supportive seating while enabling freedom to carry out their daily activities and help them get around.
- An early intervention product that focuses on simplicity, ease of use and convenience.
- Modular in construction, light to carry and simple to store.
- Supports the child in a comfortable position so they can carry out their daily activities and help with their cognitive and postural development.
- Seating system can be transferred in seconds from one base to another or remove for transportation.



## Sammons Preston Leckey® Squiggles Saddle Seat

- Designed specifically for children who just need a small element of support to help concentrate in an active seating posture.
- Saddle shaped seat is widely recognized by therapists to promote improved posture and balance in children.
- When used in combination with the pelvic harness, it will hold the child in a secure yet active position.
- Ideally suited to the clinical setting where positive positioning of the child's pelvis and feet will enable the child to focus on their therapy program.



## Sammons Preston Leckey® Advanced Seat

- Offers an unsurpassed level of positional and clinical benefits for children from a broad spectrum of varying needs.
- Works wonderfully at home, school or clinical environment.
- Myriad positional adjustments and accessory options provide a complete seating management solution.
- The height-adjustable chassis lets children fully interact with friends and family in countless activities.
- Maintains the pelvis in a neutral position with a slight anterior tilt to help in the development of a lumbar curve.
- Available in 3 sizes to accommodate children ages 1-17. Size 3 for ages 10-17 is available in Manual or Powered.



## Sammons Preston Leckey® Squiggles Stander

- Offers a versatile choice of upright, prone and supine standing all-in-one base.
- Helps children ages 6 months to 5 years old engage in wide range of activities by allowing them to convert from one standing.
- Modular in construction with two main parts - the standing support base and the chassis base.
- Standing support base can be easily clipped on and out of the chassis bases in seconds for easy storage and transportation.
- The base options range from a compact static base which folds away easily, to mobile pivot base offering a wide range of tilt options, adjusting from vertical to almost horizontal.



## Sammons Preston Leckey® Pronestander

- Offers a wide range of prone angles to allow the gradual increase of weight bearing as tolerance increases.
- Can be angled at 0, 10, 20, 30 and 40 degrees.
- Divided foot plate with sandals and Velcro foot straps are hinged for easy wheelchair transfer or emergency removal.
- Includes: Activity bowl with tray cover, flexible lateral supports, pelvic band, pelvic plate and knee supports.
- 3 sizes accommodate children through adult.



## Sammons Preston Leckey® Freestander

- This vertical stander offers multi-adjustable support and positioning at chest, hips, knees and feet.
- Design allows easy set-up during wheelchair transfer.
- Strapped sandals on the foot plate enhance custom positioning.
- Extra hip support results from the standard adjustable pelvic band and the optional pelvic plate.
- Available in 3 sizes to accommodate child through adult.
- Includes: activity bowl with tray cover, flexible lateral supports, pelvic band, knee supports, foot plate and sandals.



# Standers



## Sammons Preston Totstander

- Pre-school version of Freestander meets the needs of children at the early stages of standing.
- Accommodates children from 26" – 33" tall.



## Sammons Preston Leckey® Horizon Stander

- Versatile, attractive stander provides appropriate, gradual weight-bearing support for ages 5-18 in the prone, supine or upright position.
- Durability is featured in aluminum frame, fully adjustable supports and user-friendly powered, adjustable actuator.
- Open design offers good access from hoists and wheelchairs.
- Accessories include lateral hip and head supports and conversion kits.
- Conversion kits (Prone to Supine or Supine to Prone) provide a 90 degrees adjustment range from horizontal position to freestanding in prone and supine positions.



## Sammons Preston 3-in-1 TriStander® from Tumble Forms2

- As its name implies, the TriStander allows prone, supine or vertical standing.
- A foot-operated pneumatic tilt mechanism makes it easy to adjust from 15 to 90 degrees.
- Includes: TriStander frame, shoes with bar, hip support, headrest, knee support and trunk support.
- Each of the two sizes adjust in myriad ways to grow with the child. The first size TriStander 45 supports children 32" – 45" tall, up to 70 lbs. The TriStander 58 accommodates children 43: - 48" tall and up to 150 lbs.
- Choose from accessories such as a Push Handle, padded Adjustable Trunk Support, Head Support with Laterals, and Activity Tray. The Tri-Stander frame is also available separately.



## Otto Bock Kimba MPS Seating System

- Seat is made of incontinent-resistant foam and is covered with washable, incontinent-resistant upholstery.
- Sun canopy comes standard with an attendant window for parents or caregivers.
- Comes in two upholstery options that are child and parent friendly.
- Therapists and caregivers can make easy adjustments to the seating system while the child is still in the seat.



## Otto Bock Nurmi Neo

- Sporty design available in two fun colors: Saturn Blue Metallic and Granny Smith Green.
- Available with three different depth and height adjustable grip bars.
- Anti-tip bars provide stability to prevent tipping.
- Caster swivel wheel locks.
- Storage basket allows personal items to be kept within reach.



SAMPLE

Your Company Ad Here